

# WHAT'S NEW?

Q3 - 2008 Product Release



INTUITIVE  
SURGICAL®

Taking surgery beyond the limits of the human hand™

## EndoWrist® Graptor™ (Grasping Retractor)

New multi-functional 8mm instrument to facilitate soft tissue management

da Vinci® P/N	da Vinci® S™ P/N	Part Description	Price	Uses
400278	420278	Grasping Retractor	\$ 2,400	10 uses

Features:	Potential Benefits:
Concave fenestrated jaw profile	Allows for large, secure, atraumatic tissue purchase
Atraumatic teeth profile	Enables atraumatic grasping of bowel and other delicate tissues, as well as uterine manipulation
Precise tip grip	Designed to handle thin peritoneal tissue well
55mm jaw length with rounded edges	Facilitates atraumatic retraction of both solid and hollow viscera

### Targeted Applications (but not limited to)

- > da Vinci Sacrocolpopexy
- > da Vinci Hysterectomy for Cancer with Lymphadenectomy
- > da Vinci Low Anterior Resection



8mm EndoWrist Grasping Retractor

Required System Software: A4.3 p2 or higher

For more information,  
please call Intuitive Surgical  
Customer Service at  
1-888-409-4774

[www.intuitivesurgical.com](http://www.intuitivesurgical.com)