

Certificate

SMP GmbH
Pruefen Validieren Forschen

hereby certifies the

cleaning validation of Endowrist® Instruments

of the company

Intuitive Surgical®
950 Kifer Road
Sunnyvale, CA 94086
USA

The SMP investigation, with the project no. 05706011408, shows that the cleaning effectiveness of the Medisafe SI PCF Washer in combination with the detergent 3 E Zyme is powerful enough to assure good cleaning results of Intuitive Surgical® Endowrist® Instruments, with specific limitations:

- All instruments are visible clean
- All tips of the instruments are clean after the process.
- All tests fulfilled the acceptance criterias

Tests followed

ASTM E 2314-03
EN/ISO 15883
AAMI TIR 30



Tübingen, Germany; Oct. 30; 2006

A handwritten signature in black ink, appearing to read 'K. Roth', is written over a horizontal line.

Klaus Roth
Managing Director